



## NEW: Ready to use buffer solutions in practical Twin-neck dosing bottles

Ready to use chemsolute® buffer solutions in twin neck dosing bottles are characterized by their quality, precision and expedience in handling. The integrated calibration container allows a particularly economical dosing of solution and prevents contamination.



### Ease of use:

- Open safety lock of the integrated calibration container
- Fill calibration container with approx. 20 ml by pressing the bottle gently
- Calibrate pH electrode
- Discard used buffer solution

### Buffer solution with color coding, accuracy ( $\pm 0.01/20^\circ\text{C}$ ), traceable to NIST

	Order-No.	Unit	Price EUR/Unit
Buffer solution pH 4.00 (red)	1134.0252	250 ml	9,00
	1134.1012	1 l	20,80
Buffer solution pH 7.00 (yellow)	1137.0252	250 ml	9,00
	1137.1012	1 l	20,80
Buffer solution pH 10.00 (blue)	1140.0252	250 ml	9,00
	1140.1012	1 l	20,80



### Buffer standard solution according to DIN 19266, accuracy ( $\pm 0.010/25^\circ\text{C}$ ), traceable to NIST

	Order-No.	Unit	Price EUR/Unit
Buffer solution pH 3.776	1123.0252	250 ml	9,80
	1123.1012	1 l	21,50
Buffer solution pH 4.005	1125.0252	250 ml	9,80
	1125.1012	1 l	21,50
Buffer solution pH 7.413	1128.0252	250 ml	9,80
	1128.1012	1 l	21,50
Buffer solution pH 9.180	1131.0252	250 ml	9,80
	1131.1012	1 l	21,50
Buffer solution pH 10.012	1132.0252	250 ml	9,80
	1132.1012	1 l	21,50

# Th. Geyer

## **Deutschland**

### **Th. Geyer GmbH & Co. KG**

Dornierstraße 4  
D-71272 Renningen

Tel. +49 7159 1637-0

Fax. +49 7159 18417

E-Mail [renningen@thgeyer.de](mailto:renningen@thgeyer.de)

### **Th. Geyer Hamburg GmbH & Co. KG**

Oehleckerring 14  
D-22419 Hamburg

Tel. +49 40 533214-0

Fax. +49 40 5321004

E-Mail [hamburg@thgeyer.de](mailto:hamburg@thgeyer.de)

### **Th. Geyer - Vertriebsbüro West -**

Kirchstraße 1  
D-53797 Lohmar

Tel. +49 2246 10098-0

Fax. +49 2246 10098-29

E-Mail [lohmar@thgeyer.de](mailto:lohmar@thgeyer.de)

### **Th. Geyer Berlin GmbH**

Dillenburger Straße 53  
D-14199 Berlin

Tel. +49 30 391100-1

Fax. +49 30 3932069

E-Mail [berlin@thgeyer.de](mailto:berlin@thgeyer.de)

## **Dänemark**

### **Th. Geyer Danmark ApS**

Industrivej 6A  
DK-4000 Roskilde

Tel. +45 4630 00 30

Fax. +45 4630 00 69

E-Mail [sales@thgeyer.dk](mailto:sales@thgeyer.dk)